

Technical Solutions

AUSTRALIA PTY. LTD ACN 063 349 502

109 FERNDALE ROAD

SILVAN VIC. AUSTRALIA 3795

PHONE (03) 9737 9000 FAX (03) 9737 9111

Information Sheet

Floor Mat Alarm - Wireless

Catalogue Number: RLCB-FMA

This unit features a vinyl mat switch with a small yet powerful radio transmitter - about the size of a pack of cards. The signal is sent to a receiver buzzer box that can be located in the carer's room.

The alarm will operate over distances of up to 30 metres*.

Mat Switch

Transmitter

Receiver



Easy To Set Up:

1. Place the mat switch on the floor beside the bed.
2. Cover the mat switch with a non-slip floor mat – this will both protect and disguise the mat. **(WARNING: AN UNCOVERED MAT SWITCH CAN BE A TRIPPING HAZZARD)**
3. Fit a NEW Alkaline 9Volt battery to the Transmitter.
4. Plug the lead from the mat switch into the Transmitter and place the transmitter under the bed – out of harm's way.
5. The receiver can be moved from room to room to suit the needs of the carer/s. More than one receiver can be set-up to respond to the call if required.
6. Simply plug the receiver power supply into a 240Volt power point.
7. The receiver has two modes of operation:
 - **Instant** - The buzzer will sound only while someone is standing on the mat switch.
 - **Latch** - The buzzer will start when someone stands on the mat switch and will keep buzzing till the carer resets it.

IMPORTANT Maintenance:

- Test the transmitter battery regularly with a battery tester (available from Technical Solutions) Replace when necessary, with a NEW Alkaline 9Volt battery. Battery life will vary with frequency of use. If radio range starts to decrease, it is likely that the battery is going flat.
- Take care of your mat switch. Do not fold, as this will destroy the switch. Store flat do not roll. Avoid placing sharp objects on mat switch.

* The range will be dependant on environmental factors.